

**MEDICUBE DIGITAL AGE-R**



**AGE-R DEVICE  
BOOSTER-H**



**USER  
MANUAL**

## Precautions During Use

1. Dropping the product during use or lifting the product by holding only the cap may cause damage and product failure.
2. Make sure no liquid or cosmetic product enters the device.
3. Contact the manufacturer if the product does not work; do not disassemble it.
4. Make sure that the device and the charging cable are connected properly.
5. Before use, make sure that the cables and any accessories are connected properly.
6. Remove any metal accessories that may come into contact with the user's body.
7. Do not use together with electrical appliances such as electric mats.
8. Do not use next to a fire or in hot places, such as under direct sunlight. Keep body and cable away from heat sources.
9. Do not hold the device in one area for long periods of time. Continuously massage in motions.
10. Beware of contact with other products as the device may malfunction.
11. To prevent overheating of internal components, turn off the device for 10 minutes when using the product more than 3 times in a row. Continuous use may shorten life of the device and may cause malfunctions.
12. Do not use it for any purpose other than its intended purpose. It may be the cause of accidents.
13. Please do not use on persons who are paralyzed or unconscious.
14. Do not tap the roller head on solid surfaces, as this may cause malfunction.
15. Do not hit the device against hard objects or drop the device on the floor.
16. Do not use in high temperatures, high humidity, or in areas with excessive dust, such as bathrooms, kitchens, and cars.
17. When using, stay clear of electrical appliances such as other beauty equipment, TV, radio, etc.
18. Do not use the device around the Adam's apple.
19. Do not use on any voltage level other than the voltage level it is rated for. Using above the recommended level may cause excessive heat, fire, short circuit, or malfunction.
20. Do not handle the device and charging cable with wet hands (this may cause electric shock).
21. Do not repair, disassemble or modify the device on your own (this may cause product damage).
22. The device does not operate when charging.

## Precautions for User

Please do not use the device if you are in these condition ; particularly bad health condition, drunken condition, taking sleeping pills, feeling severe fatigue and being treated by doctor. If you feel pain or health problem during use, please stop using. If you use the device on same area repeatedly or use in wrong method, it may damage the skin. If you feel something wrong with your skin, keep the affected area cold and consult with specialist. \* If your skin is sensitive, the skin can get slightly red even in general use. Please well noted with the manual and precaution and go through appropriate after care.

### Do not use if you are in condition down below

1. Anyone whose skin is red and tingling due to excessive tanning.
2. Anyone whose skin is specially sensitive due to allergic constitution.
3. Anyone who is using an electronic device directly worn in body such as cardiac pacemaker.
4. Anyone who is using an electronic medical device for life retaining such as artificial heart-lung machine.
5. Anyone who is using a mounted electronic medical device such as ECG device.
6. Anyone who is using a steroid type hormone drug for long time.
7. Anyone who has telangiectasia due to liver function disorder.
8. Anyone who is outpatient receiving treatment.
9. Area where acne festering and inflamed.
10. Dermatitis area caused by cosmetic product, anyone who has atopic dermatitis.
11. Freckle area accompanying itchiness, freckle caused by physical stimulation.
12. A part of the body that contains metals, plastics, silicon, etc.
13. Hemophilia, heart disease, acute disease, hemorrhagic area, malignant tumor, fracture area, infectious disease, eyelid area, skull area. Children, pregnant women, lactating women and elderly.

\* Since the device is not a medical device, do not use for therapeutic purposes.

## Precaution During Storage

1. Store at place where it is hard for children to reach and where there is no moisture.
2. Do not store at place with magnetic or electromagnetic wave occurrence. It may cause breakdown.
3. After using the head part, wipe it with wet tissue or clean gauze, soft cloth, etc. And store it in dry condition.
4. Store the body part with the power turned off.
5. Do not wipe the product with inflammable material such as thinner and benzene. It can cause discolor and deformation due to a chemical reaction.

## Specification

Model Name	ME-R5TH	Product Name	MEDCUBE AGE-R BOOSTER-H
Manufacturer	Easystem Co., Ltd.	Contact	031-312-1305
Address	(Prime Castle Knowledge Industry Center) #304, 25, Podowowo 116beon-gil, Siheung-si, Gyeonggi-do, Republic of Korea		
CE Certificate Registration No.	R-R-pg1-ME-R5TH	Method of Use and Precautions	Refer to the user manual
Manufacturing No., Manufactured Date	Refer to the specification label	Other Necessary Matters	Refer to the user manual

### Product Specification

#### Product Dimensions:

158.2\*149.7\*43.0mm / 168.2\*159.7\*36.0mm

Input Voltage 5V/1A (5.00A)

Battery

450mAh Lithium-polymer

Power Consumption

maximum of 1W

Weight 75g (±1.7%)

Charge Time about 3 hours

MADE IN KOREA

### Charger Specification

Quality Indication in accordance with the Electrical Appliances And Consumer Products Safety Control Act Safety Confirmation Report No. XJ100743-100010



Model Name DTP-100030

Product Name Battery

Rated Voltage 5.7V, 450mAh

Manufactured Date: Indicated separately at the bottom of the product

This device observes Article 15 of the Federal Communications Commission (FCC) regulation and fulfills two conditions below.

1) Non-creation of harmful substances designated by the relevant laws.

2) Absorbing power of exterior impact that may cause malfunction not intended by the user.



### Exchange / Return Policy

#### 1. Warranty

For all Age-R devices, we offer a 1-year limited warranty period from the date of your original purchase. If your device does not work properly due to a manufacturer's defect, we will replace the defective device free of charge within 1 year. This warranty does not cover the device damaged by: 1) Accident, misuse, or alteration 2) Use with unauthorized accessories or use other than instructed, 3) Connecting to incorrect current and voltage (We only provide C-type USB cable) This warranty only extends if the device are used in conjunction with authorized device accessories (such as the enclosed charging cables).

#### 2. Returns

For the Age-R device, all returns must be initiated within 30 days of receiving your order. Please refer to the information as shown below. Simple Return: A refund due to a change of mind is only applicable within 30 days of receiving and will only be accepted if the product is returned in original, sealed, and unopened condition. In the case of a return for simple reason (change of mind), the customer needs to bear the shipping fee. For international orders: You will be responsible for paying for customs fee and returning shipping fee. Additionally, we cannot provide free return shipping service for international customers.

Damaged or defective product: 1) Please reach out to us via our Customer Service Center (Email with the following details: First name and last name, email address, phone number, order number. Kindly attach a photo and video describing your issue. 2) We will proceed to issue a confirmation email for you (Please note that our Customer Service usually takes 2-3 business days to respond). 3) After receiving a confirmation email, please send the product back to our warehouse via the return address. You will be informed of the return address after the email confirmation. 4) Once the product has arrived at our warehouse, our staff will inspect the condition of the product. 5) In case of product defect, we will proceed to exchange and send a new product to you. \*Note: This process may take up to 1 month, including the time taken to ship the product back to our warehouse, and the time taken to inspect the product.

### Subject of Recommendation of Lithium Battery Discarding and Charging Method

Warning : Since it may accompany danger such as heating, fire, and explosion, please observe the following matters. 1) Handling Matters Since a swollen battery that can be visually identified may be dangerous, immediately contact the manufacturer or the seller. Use only the designated charger. Keep away from fire (Do not put inside a microwave oven). Do not leave it inside a car in summertime. Do not store or use at high temperature and humid place such as saunas. Do not use it for a long time while placing it on blanket, electric pad, and carpet. Do not store at a sealed space while leaving the power on. Be cautious not to contact the electrical terminal with metal products such as necklaces, coins, keys, and watches. Use the secondary lithium battery warranted by the manufacturer. Do not implement behaviors such as disassembly, compressing, and penetrating. Do not apply abnormal impact such as dropping it from a high place. Do not expose to a temperature above 50 Celsius. Do not contact with humidity. 2) Disposal Guideline Do not arbitrarily dispose used up battery. When disposing the battery, do not incinerate or put it inside a fire. Explosion and fire may occur due to heat. The method of battery disposal may vary according to countries and regions. Dispose with appropriate method. 3) Recommendations regarding Charging Method Do not charge the battery with the method other than the one shown in the user manual. Fire or explosion may occur. Complete Charge Time : About 3 hours \*Charge time may vary depending on the battery remnant amount of portable battery or surrounding temperature or charging environment. It may change without notice for performance improvement in specifications and exterior appearance.

## Description

---



## Operation Method

---

1. Turn on the product by pressing ③ Power/Level Button for 2~3 seconds.  
\* When turning on the power, ② Color LED Display blinks with a sound "ON".
2. The product is automatically set to STEP 1 with a voice alarm after turning on. Adjust the level by pressing ③ Power/Level Button. Each time the button is pressed, the level goes up by one level and ② Color LED Display blinks with a level setting voice. \*It is able to set from Level 1 to Level 5. Please make sure to use it from Level 1 for initial use.
3. Once the level setting is complete, make sure that the front part of the head is fully attached to the skin. \*When you use it on narrow area such as next to a nose, smile lines, use the protruding part on the upper part of the head. \*For [Low to Boost], please refer to the details in the User Guide.
4. The operating time will be informed by the voice alarm every 1 minute and the product will be automatically turned off when it comes to 5 minutes.  
\*If you want to stop operating, press ③ Power/Level Button for 2 seconds or more. The product will be turned off with a sound "OFF".  
\*If the battery is low, ② Color LED Display will blink in red and the product will be automatically turned off with a voice alarm of "Low Battery".
5. Clean ① Head with a dry cloth or cotton pad after using for hygienic use.
6. Connect the Type C charger enclosed in the product box to the bottom of the rear to charge the device.  
\* ② Color LED Display turns on in red while charging and in green when fully charged.  
\* The product does not work while charging.

## Contents

---



**Medicube AGE-R Booster-H**  
USB Charging Cable (USB-C)  
User Manual

**DEDICATED TO  
REAL RESULTS  
MEDICUBE**



**STEP 01**

After washing your face gently, applying the skincare product on your face.



**STEP 02**

Fully attached the head to the the skin surface and rub the skincare product to be absorbed into the skin.



**STEP 03**

Refer to the boosting method, please use the device until the skincare product on your face is fully absorbed.



**STEP 04**

When the skincare product is fully absorbed, finish with intensive care.

## Bring out the natural glow of the skin with AGE-R's deep absorption technology



After cleansing, apply skincare product thoroughly. Massage gently with the head of device for boosting effect **in closed contact** with the skin surface.

### Boosting Method

---



#### Watering

Push the head of device against the skin surface and rolling it like a drawing the circle on the skin.



#### Boosting

Push the head of device against the skin surface and sliding it on the skin.



#### Focus

Push the head of device against the skin surface and pressing it one spot at a time on the skin.

### Read before use!

---

It is recommended to start from level 1 when using it for the first time. Please adjust the level appropriately while continuously using it.



#### Caution

Please avoid using uvula (Adam's apple), wounds, and surgical area.



#### Recommended starting level

Start slowly from the lowest setting and find the appropriate level while continuously using it.



#### BOOSTER-H Use Cycle

You may use the Ago-R Booster-H 1-3times a day. Please gradually increase the frequency of using.



## **MEDICUBE DIGITAL AGE-R**

**So that everyone can enjoy the digital service  
which was privilege for few people**



**AGE-R DEVICE  
BOOSTER-H**

**medicube** 