

MEDICUBE CLINIC AT HOME



AGE-R DEVICE
DERMA EMS SHOT



USER
MANUAL

Description



Operating Guide

1. Press and hold the ⑥ Power button to turn on the power.
*When the device is turned on, the remaining battery level will be displayed on the ③ LED Indicator in white.
2. Press and hold the ④ Mode/Intensity adjustment button to change the mode.
● Skin ● Up ● Body
3. After selecting the mode, adjust the intensity level by pressing the ④ Mode/Intensity adjustment button quickly. When the ③ LED indicator is on level 1, one LED light is lit. An additional LED light is then lit every time the level goes up. (The device can be set between levels 1-4)
4. Apply adequate amount of gel to the desired area and gently massage the ① Roller Heads against the skin in an upward motion.
*The notification rings every 5 seconds when the device is in operation.
5. The device turns off automatically after 10 minutes of usage. To turn off the product during use, press and hold the ⑥ Power button.
6. After use, take a soft cloth or cotton pad to clean the ① Roller Heads before storing it.
7. Charge the device using the C-type charger enclosed in the box. When charging, the battery level will be displayed on the ③ LED Indicator in white.
*The device cannot be operated while charging.



Contents

- Medicube Age-R Device
- USB Charging Cable
- User Manual

Let us focus on your skin,
so you can focus on yourself.

Cautions

1. Dropping the product during use may cause damage and product failure.
2. Make sure no liquid or cosmetic product enters the device.
3. Contact the manufacturer if the product does not work; do not disassemble it.
4. Make sure that the device and the charging cable are connected properly.
5. Before use, make sure that the cables and any accessories are connected properly.
6. Remove any metal accessories that may come into contact with the user's body.
7. Do not use together with electrical appliances such as electric mats.
8. Do not use next to a fire or in hot places, such as under direct sunlight. Keep body and cable away from heat sources.
9. Do not hold the device in one area for long periods of time. Continuously massage in upward and downward motions.
10. Beware of contact with other products as the device may malfunction.
11. To prevent overheating of internal components, turn off the device for 10 minutes before reusing. Continuous use may shorten life of the device, may cause accidents.
13. Please do not use on persons who are paralyzed or unconscious.
14. Do not tap the roller head on solid surfaces, as this may cause accidents.
15. Do not hit the device against hard objects or drop the device on the floor.
16. Do not use in high temperatures, high humidity, or in areas with excessive dust, such as bathrooms, kitchens, and cars.
17. When using, stay clear of electrical appliances such as other beauty equipment, TV, radio, etc.
18. Do not use the device around the eyes or around the Adam's apple.
19. Do not use on any voltage level other than the voltage level it is rated for. Using above the recommended level may cause excessive heat, fire, short circuit, or malfunction.
20. Do not handle the device and charging cable with wet hands (this may cause electric shock).
21. Do not repair, disassemble or modify the device on your own (this may cause product damage).
22. The device does not operate when charging.

Warnings

Do not use if you are in poor health, under the influence of alcohol, taking sleeping pills, suffer from chronic fatigue, or are currently being treated by a doctor for a serious medical condition. Stop use immediately if you experience bodily pain or abnormalities during use. If used incorrectly or repeatedly on the same area, it may damage your skin. If any skin abnormalities occur, immediately cool the area and consult a specialist.

*For those with sensitive skin, the skin may turn slightly red even during normal use. Read the manual carefully and perform proper aftercare.

Do not use Derma EMS Shot Age-R Device if you:

1. Are currently suffering from sunburn or have red, inflamed skin.
2. Have allergies or especially sensitive skin.
3. Have a pacemaker, or any other active electrical implant anywhere in the body.
4. Are using electronic medical devices for life support, such as a heart/lung machine.
5. Are wearing electronic medical devices, such as an ECG device or a heart monitor.
6. Are a long-term user of steroid hormone drugs.
7. Have impaired liver function with capillary abnormalities.
8. Are undergoing outpatient treatment.
9. Have inflamed acne or boils.
10. Have dermatitis due to cosmetic application, or atopic dermatitis.
11. Have melasma accompanied by itching or melasma caused by physical irritation.
12. Have metals, plastics, silicon etc. implanted in the body.
13. Have hemophilia, heart disease, acute disease, hemorrhage, malignant tumor, fracture, infectious disease, or physical trauma. Avoid using on the eyelid and skull areas.

Maintenance

1. Keep out of reach of children and store the product in a dry place.
2. Do not store the product near magnetic or electromagnetic waves. This may cause damage to the product.
3. After use, wipe the roller head part with wet wipes, clean gauze or soft cloth, and keep it dry.
4. Be sure to turn off the power before storing.
5. Do not use flammable substances such as thinner or benzene to clean the product, as this may lead to chemical reactions, and cause discoloration and damage.

Specification

Equipment Name	MEDICUBE AGE-R DEVICE (Model Name: ME-AGE-R)	Model Number	ME-AGE-R-01
 Manufacturer	BSL CO., LTD	Contact	070-4030-6810
		Manufacturer Address	301, 302, 9-21, Dongtansandan 9-gil, Hwaseong-si, Gyeonggi-do, Republic of Korea
KC Certificate Number	R-R-pg1-ME-AGE-R-01	Serial Number, Manufacturer Date	Refer to the product label
Directions and Precautions	Refer to the manual	Other necessary information	Refer to the manual
Product Specifications Product Dimensions: 172*70*45mm Input Voltage: DC5V / 2A Battery Capacity: 3.7 V / 500 mAh Power Frequency: 1KHZ Weight: 240 g Battery Charging Time: Approx. 2 hours		Battery Specifications Quality Indication in accordance to the Electrical Appliances and Consumer Products Safety Control Act Registration No.: ZU10701-17002J Model Name: YJ002535 Product Name: Lithium Polymer Battery Battery Rating: 3.7 V, 500 mAh Manufacturer: YJ POWER GROUP LIMITED  Manufacturing Country: China Manufacturing Date: Marked on the bottom of the product	
MADE IN KOREA			

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.



R-R-pg1-ME-AGE-R-01

Exchange / Return Policy

1. In case of product defect, Medicube AGE-R Device™ can be exchanged within 30 days of received date. The exchange process is outlined as below. 1) Send us an email to the cs_email on the website you purchased with the following details: First name and last name, email address, phone number, order number. Attach a photo or video describing the issue. 2) We will proceed to issue a confirmation email (our Customer Service usually takes 2-3 business days to respond). 3) After receiving a confirmation email, send the product back to our warehouse via the Return address. Return address will be shared after the email confirmation. 4) Once the product has arrived at our warehouse, our staff will inspect the condition of the product. 5) In case of product defect, we will proceed to exchange and send a new product to the customer. *Note: This process may take up to 1 month, including the time taken to ship the product back to our warehouse, and the time taken to inspect the product. 2. Medicube will not refund any products returned due to a "change of mind" or "simple remorse". If returned products do not fulfill the requirements for exchange, they may not be refunded and Medicube will notify the customer if this is the case. The customer will be asked if they want the products to be sent back to them, but the customer will be responsible for any delivery costs. If the customer does not pay for the delivery back within 30 days, Medicube will deem the products to have been abandoned by the customer, including any right of refund.

This warranty will not be re-issued, please retain this copy.

[Instructions for lithium battery disposal & recharging]

WARNING: Please read the following instructions to avoid risk of overheating, fire, explosion, etc. 1) **How to handle** Contact the manufacturer or supplier immediately if batteries appear visibly swollen as it is a dangerous hazard. Keep away from fire or flame (do not put in the microwave). During summer months, do not leave unattended in a car. Do not use or store in hot, humid areas such as saunas, etc. Do not use for extended periods of time while placed over blankets, electric mats, carpet, etc. Do not store for extended periods of time in enclosed spaces while the power is turned on. Be careful not to touch the battery terminals with metal objects such as needles, coins, keys and watches. Use only secondary batteries verified by the manufacturer. Do not disassemble, puncture or compress. Avoid dropping from high places or other forms of impact. Do not expose to temperatures above 50°C. Avoid exposure to humidity. 2) **Battery Disposal** Do not throw away used batteries. Do not incinerate or place in fire to dispose. High heat may cause fire or explosions to occur. Battery disposal methods vary by country and region. Please follow appropriate local battery disposal methods. 3) **Instructions for recharging** Follow the instructions to recharge. Fire or explosion may occur if the product is recharged with other methods not stated in the instructions. Time for full charge: Approx. 2 hours *Charging time may vary according to the remaining battery level, ambient temperature, or charging environment of the portable battery. Specifications and appearance are subject to change without prior notice to improve performance.

Quick Start Guide

Face Slimming → **SLIM**

Face Up → **UP**

Body Contouring → **BODY**



STEP 01

After cleansing and applying toner, apply the Booster Gel to the roller heads of the Age-R Device.

***Gel must be applied in order to feel the energy.**



STEP 02

Select desired mode and intensity level. ***When using the product for the first time, start at the lowest setting.**

***When the device is in operation, it will emit a notification sound every 5 seconds.**



STEP 03

Apply the device to the desired area and massage in an upward direction.

***Recommended order of use: Jawline → Laugh Lines → Cheeks → Crow's Feet → Temples → Forehead**



STEP 04

After 10 minutes of use, pat the skin to allow the gel to fully absorb.

***If the skin feels dry, apply a moisturizer to finish.**

Getting the Best Results from Your Age-R Device



Jowl & Jawline

Massage from the tips of the chin to the bottom of the ear, pulling the skin gently in an upward motion.



Smile Lines

Divide the area into 2 horizontal zones, and massage both areas in an upward motion.



Cheek Care

Divide the area into 3 vertical zones, and massage each area in an upward motion from the jawline to just below the eye.



Crow's Feet

Gently roll in an upward motion to relax the muscles.



Temples & Forehead

Divide the area into 3 vertical zones, and massage each area in an upward motion.



Shoulder & Trapezius Muscle

Massage the muscles by rolling back and forth in upward and downward motions.

Mode Settings

Each mode designed to replicate the feeling of an esthetician's acupressure massage techniques.



Slim Mode



Up Mode



Body Mode

Mode Description

Mode

SLIM

Slim

Stimulates face and helps with slimming the face



*This mode is designed to simulate the Tapotement massage technique.

UP

Up

Help to stimulate face to improve the appearance of the face's contours



*This mode is designed to simulate the Petrissage massage technique

BODY

Body

Helps to relax and improve the appearance of the body's contours



*This mode is designed to simulate the Friction massage technique.

Read Before Use

The Age-R Device works by moving and massaging your skin.



Caution

Avoid using around the eyes or around the Adam's apple.



Caution

Avoid using on the armpits and the chest area.



Test Before Use

Before use, apply the gel and test the device on the palm of your hand (by placing both the roller heads on the spots as shown above). Start the device from the lowest setting.

- When using the device for the first time, start slowly from the lowest setting.
- It is highly recommended to test the device before using on the face.

*Both roller heads must maintain contact with the skin simultaneously for the medium-frequency-energy to be emitted successfully.

Medicube presents,
A new good energy beauty massaging device for better skin care condition.

We promise the highest quality for
your safe, daily use of the Age-R Device.

Certified by safety federations around the world (Europe, USA, Japan and Korea)

The Medicube Age-R Device has been approved by the following 4 global certification organizations: KC, CE, FCC, and PSE. It is also registered with the US FDA.



QUALITY YOU CAN TRUST

Designed with medical-grade surgical steel, the very same material you would find in a surgery room, you can trust it to be safe, durable, and resilient to rust.

MADE IN KOREA

The Medicube Age-R Device is designed, manufactured and assembled in Korea, a global leader in cosmetics and beauty innovation.



FAQ

Q1.
Do I really need to use the Booster Gel with the device?

The Age-R Device is a beauty massaging device that operates by delivering energy, and therefore must be used with products that are capable of delivering the energy. It is highly recommended that you use the dedicated gel serum that was developed alongside the device for best results. Medicube's Booster Gel contains plant stem cell extract, which helps the improvement of skin condition.

Q2.
How much gel do I need to apply?

We recommend applying a generous layer on the desired area, as well as a small dollop on the roller heads of the device.

TIP: Using more gel strengthens the sensation of stimulation. If you feel that the strength of the device is too weak or too strong, try adjusting the amount of gel applied.

Q3.
I feel my teeth aching when I use it around my nose and mouth area, is that normal?

Due to the deleving energy by massaging around the mouth region, you may feel a chilly and tingling sensation in your teeth. You may also feel this on your scalp when you use the device around your forehead. These symptoms are normal and commonly occur even when receiving treatments at aesthetic shops as well. The tingling sensation is temporary, and usually disappears within approximately 10 minutes.

Q4.
Should I wash off the gel once I'm finished?

Our gel was developed to be multifunctional as both a current conductor as well as a serum that helps to improve the facial contours through moisturization of the skin. You can wipe it off if you would like, but we recommend patting it into the skin after finishing your massage. If the skin feels dry even after applying the gel, apply your favorite moisturizer on top to finish.

Q5.
Can it be used on areas other than the face?

It can be used on various parts of the body, such as the shoulders, calves, and upper arms. Do not use around the eye area, Adam's apple, skull, armpit, and chest regions. Do not use on areas with cuts or open wounds of any type.

Q6.
Do I really have to use it for 10 minutes every day?

Due to the nature of the device, using it for an extended period of time in one sitting can cause irritation. This is why we recommend using it within the time limit of 10 minutes per day. **TIP:** There may be individual differences in the intensity of sensation. Please adjust the usage time according to your comfort level.

Q7.
How can I clean the device?

After use, wipe down the roller heads with a wet wipe or cosmetic cotton pads. Store it in a dry area.

Q8.
The device won't turn on at all!

If the battery is insufficient, the device will be unable to operate. In this case, charge the device sufficiently and resume use.

Caution

- Avoid contact with the eye area, Adam's apple, armpit, and chest regions.
- If you have received laser treatments, chemical peels, or plastic surgery of any kind, consult with your doctor before using the Age-R Device.
- Before use, it is recommend to remove all metal accessories and jewelry.
- Take caution when handling the device, as it can be easy to drop due to the weight of the device.